

TO,


Subject :- Submission of Quotation for the Supply of Electronic Tourniquet
to Dean Govt. Medical College, Akola .

Dear Sir,

You are requested to submit your scaled quotation for the supply of
Machine/Surgical Item to this office with the terms and conditions mentioned
below

- 1)The quotation should be send in scaled cover (duly sealed wax) by the post
or by hand delivery so as to reach this office on be / / 2018
- 2) Quotation will be opend in the office of the dean on / / 2018
- 3)Quotation which are received late will not be accepted and any circumstances.
- 4) The words quotation for and last date , must be mentioned on envelope .
- 5) Rate should be for free delivery at the Dean , Govt. Medical college , Akola.
premises only .
- 6)Quoted good should be strictly according to the make & the specificationasked in
enquiry, and make of item by you should be stated in your quotation .
- 7)In case the order are placed with you , the order will have to be executed in full
within the stipulated time . If you fail to comply , the orders may be cancelled
and you will be declared defaulter , no enquiry will be sent in future without
assigning any reason.
- 8) The undersigned reserves the right to accept or reject any or all quotation. . The
quotation should be sent in the name of the / Dean , medical college & Hospital,
Akola.
- 9) Please quote the sale tax /GST No. registration no . in your quotation

Sr.no	Name of Reagent	Qty
1	Electronics Tourniquet	
	(As per Specification)	


Dean
Govt. Medical College,
Akola

Specification of Tourniquet Machine

Cuff pressure range 10 to 430 mmHg.

Pressure regulation +10mm hg of set point

Online setting Increase and decrease in pressure

Timer Range from 9hrs to 59 minutes,
Decrement mode form set value.

Intemal least count One minute, Internal 1/1000 second

Memory Pressure set in earlier session is
Stored and displayed, when machine
Is swiched ON again,

Pressure Regulation Maintains set Presure under leak conditions.

Backup 3 Hour battery backup (on full charge)
Switchover Automatic switchover from mains to
Battery mode.

Power 230v (AC)/50Hz +10%
Stabilizer Inbuilt

Digital display Digital display of time elapsed, set time,
Actual and set Pressure.

Cuffs Different sizes of 5 cuffs-washable & easy
Fitting (pediatric, small,Medium, Big & large.

Flush instantenous Inflate & Deflate.

Body Plastic

Alarm Activation along with display Indication for the following:

1. After set time is over it will show t.Ou,
2. Only once it will show P.HI, when pressure goes 25mm above set pressure
3. Only once it will show P.Lo, when pressure goes 25mm below set pressure
4. In battery mode, when battery voltage reaches 9V, it will show bat, Lo.